



Electrosurgical Unit For VET




● PPS (Patient Protection System)

PPS applied to monopolar with micro pulse technology makes the returnpad burn eliminate by 500 times/sec contact monitoring system between patient and returnpad. If the contact condition is not good enough to safety level, PPS controls the output power automatically.

● ACW (Arc-Controlled Waveform)

New technology ACW has absolute output power at every mode of Coagulation, Blend, Bipolar like variable power of general ESU. However unlike general ESU which has "off time" on waveform, it has variable frequencies for pulse continuity. The feature keeps the output power of general ESU strong and makes the electric arc be dramatically depleted. Finally it minimizes not only the carbonization but also second bleeding by adhesion of the tissue.

● ACW solution



Specifications

	HF - 200	HF - 300
Electrical Specification		
Supply voltage	220V	220V
Frequency	60Hz	60Hz
Power input during max.	900W	1100VA
Fuse	T10A	T10A
Power Output Specification		
Maximum Cut[555KHz]	200W	300W
Maximum Blend[ACW]	150W	200W
Maximum Coag[ACW]	150W	200W
Maximum Bipolar[ACW]	70W	100W
Special Function		
Safety system	PPS(Patient Protection System)	
New Technology	ACW(Arc Controlled Waveform)	
Demensions and Weight		
Width x Depth x Height	330 x 350 x 130	330 x 350 x 130
Weight	7.5Kg	7.5Kg

經銷商

悅田有限公司

234 台灣新北市永和區中正路 718 號 8 樓

電話:(886)286603155 傳真:(886)229279544

電郵:witech@ms71.hinet.net

thomas@witech.com.tw

網址:witech@ms71.hinet.net